




GOAT ANTI-TLR6 ANTIBODY





SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Aliquot and store at -20°C. Minimize freezing and thawing. Instructions

Synonym /

TLR6|toll-like receptor 6|OTTHUMP00000158730|CD286 Alias

Names

Immunofluorescence: Strong expression of the protein seen in THP-1 and A549 cells. Usage Recommended concentration: 10µg/ml. Flow Cytometry: Flow cytometric analysis of **Summary**

THP-1 cells. Recommended concentration: 10ug/ml.

Accession

ID


NP_006059.2

Blocking Peptide

EBP10208









Telephone: (650) 697-3600

Peptide with sequence C-RGSTLLNFTLNH, from the internal region of the protein sequence according to **Immunogen**

NP_006059.2.

Peptide Sequence

C-RGSTLLNFTLNH

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping Instructions

Refrigerated

Predicted

Human, Cow, Pig, Sheep **Species**

Reactive

Species

Human

Human

10333

Gene ID

Product Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

ELISA

Detection

Antibody detection limit dilution 1:1000.

Limit

Application

Pep-ELISA, FC, IF **Type**

GALLERY IMAGES

