




# GOAT ANTI-TLR6 ANTIBODY

SKU: EB10208



## SPECIFICATIONS

**Formulation** Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

**Unit Size** 100 µg

**Storage Instructions** Aliquot and store at -20°C. Minimize freezing and thawing.

**Synonym /**

**Alias Names** TLR6|toll-like receptor 6|OTTHUMP00000158730|CD286

**Usage Summary** **Immunofluorescence:** Strong expression of the protein seen in THP-1 and A549 cells. Recommended concentration: 10µg/ml. **Flow Cytometry:** Flow cytometric analysis of THP-1 cells. Recommended concentration: 10ug/ml.

**Accession ID** NP\_006059.2

**Blocking Peptide** EBP10208



<b>Immunogen</b>	Peptide with sequence C-RGSTLLNFTLNH, from the internal region of the protein sequence according to NP_006059.2.
<b>Peptide Sequence</b>	C-RGSTLLNFTLNH
<b>Purification Method</b>	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
<b>Shipping Instructions</b>	Refrigerated
<b>Predicted Species</b>	Human, Cow, Pig, Sheep
<b>Reactive Species</b>	Human
<b>Human Gene ID</b>	10333
<b>Product Grade</b>	<a href="https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png">https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png</a>
<b>ELISA Detection Limit</b>	Antibody detection limit dilution 1:1000.
<b>Application Type</b>	Pep-ELISA, FC, IF

## GALLERY IMAGES

