


# GOAT ANTI-TLR7 ANTIBODY

**SKU:** EB07806



## SPECIFICATIONS

**Formulation** Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

**Unit Size** 100 µg

**Storage Instructions** Aliquot and store at -20°C. Minimize freezing and thawing.

**Synonym /**

**Alias** toll-like receptor 7|TLR7

**Names**

**Accession ID** NP\_057646.1

**Blocking Peptide** EBP07806

**Immunogen** Peptide with sequence C-DNHVAYSQVFKET, from the C Terminus of the protein sequence according to NP\_057646.1 .

**Peptide Sequence** C-DNHVAYSQVFKET

**Purification Method** Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

|                              |                                                                                                                                                                                                                                   |
|------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Shipping Instructions</b> | Refrigerated                                                                                                                                                                                                                      |
| <b>Predicted Species</b>     | Human, Mouse, Rat, Dog                                                                                                                                                                                                            |
| <b>Reactive Species</b>      | Human                                                                                                                                                                                                                             |
| <b>Human Gene ID</b>         | 51284                                                                                                                                                                                                                             |
| <b>Product Grade</b>         | <a href="https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png">https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png</a> |
| <b>IHC Results</b>           | Paraffin embedded Human Tonsil. Recommended concentration: 5µg/ml.                                                                                                                                                                |
| <b>ELISA</b>                 |                                                                                                                                                                                                                                   |
| <b>Detection Limit</b>       | Antibody detection limit dilution 1:64000.                                                                                                                                                                                        |
| <b>Application Type</b>      | Pep-ELISA, IHC                                                                                                                                                                                                                    |

## GALLERY IMAGES

