


Telephone: (650) 697-3600

GOAT ANTI-TLR7 ANTIBODY

SKU: EB07806



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 μα

Storage

Aliquot and store at -20°C. Minimize freezing and thawing. Instructions

Synonym /

toll-like receptor 7|TLR7 Alias

Names

ID

Accession

NP_057646.1

Blocking

Peptide

EBP07806

Peptide with sequence C-DNHVAYSQVFKET, from the C Terminus of the protein sequence according to **Immunogen**

NP_057646.1.

Peptide

C-DNHVAYSQVFKET

Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.





Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

Shipping Refrigerated Instructions

Predicted

Human, Mouse, Rat, Dog

Species Reactive

Human **Species**

Human **Gene ID**

51284

Product

 $https://prod-vector-labs-pimcore-assets.s 3. us-east-1. a mazon aws.com/assets/products/image/elite_medium.png$ Grade

IHC Results Paraffin embedded Human Tonsil. Recommended concentration: 5µg/ml.

ELISA

Antibody detection limit dilution 1:64000. Detection

Limit

Application Pep-ELISA, IHC Type

GALLERY IMAGES

