



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600


GOAT ANTI-TMPRSS2 ANTIBODY

SKU: EB08219



Telephone: (650) 697-3600







Telephone: (650) 697-3600

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 μg

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias

TMPS2|epitheliasin|PRSS10|transmembrane protease, serine 2

Names

Accession

NP_005647.2

Blocking

EBP08219 **Peptide**

Peptide with sequence C-TDWIYRQMRADG, from the internal region of the protein sequence according to **Immunogen**

NP 005647.2.

Peptide Sequence

C-TDWIYRQMRADG

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Human **Species**

Reactive

Human **Species**

Human

7113 Gene ID

Product

 $https://prod-vector-labs-pimcore-assets.s 3. us-east-1. a mazon aws.com/assets/products/image/elite_medium.png$ Grade

IHC Results Paraffin embedded Human Kidney. Recommended concentration: 6-8µg/ml.

ELISA

Detection Antibody detection limit dilution 1:8000.

Limit

Approx 55-60kDa band observed in Human Pancreas and Kidney lysates (calculated MW of 53.9kDa according Western

to NP_005647.2). This molecular weight is routinely observed by other sources. Recommended concentration: Blot

0.1- $0.5\mu g/ml$. Primary incubation 1 hour at room temperature.

Application

Pep-ELISA, WB, IHC Type

GALLERY IMAGES





Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

