

GOAT ANTI-TNNI1 (AA82-96) ANTIBODY

SKU: EB12039



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias Names	troponin I, slow-twitch isoform troponin I, slow skeletal muscle TNN1 SSTNI troponin I type 1 (skeletal, slow) TNNI1
Accession ID	NP_003272.3
Blocking Peptide	EBP12039
Immunogen	Peptide with sequence DIEAKCLHNTREIKD, from the internal region of the protein sequence according to NP_003272.3.
Peptide Sequence	DIEAKCLHNTREIKD
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Dog, Pig, Cow
Reactive Species	Human, Rat
Human Gene ID	7135
Mouse Gene ID	21952
Rat Gene ID	29388
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:8000.
Western Blot	Approx 25kDa band observed in Human Skeletal Muscle lysates and approx. 24kDa in Rat Skeletal Muscle lysates (calculated MW of 21.7kDa according to Human NP_003272.3 and Rat NP_058880.1) Recommended concentration: 0.2-0.6µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

