


GOAT ANTI-TRANSCRIPTION FACTOR E2F1 (AA314-327) ANTIBODY

SKU: EB12261



250kDa
150kDa
100kDa

75kDa
50kDa
37kDa

25kDa
20kDa
15kDa

SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	transcription factor E2F1 retinoblastoma-binding protein 3 retinoblastoma-associated protein 1 RBBP-3 RBAP-1 PRB-binding protein E2F-1 PBR3 RBP3 RBBP3 RBAP1 E2F-1 E2F transcription factor 1 E2F1
Accession ID	NP_005216.1
Blocking Peptide	EBP12261
Immunogen	Peptide with sequence C-QEVTSEENRATDS, from the internal region of the protein sequence according to NP_005216.1.
Peptide Sequence	C-QEVTSEENRATDS
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	1869
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA	
Detection Limit	Antibody detection limit dilution 1:128000.
Western Blot	Approx 60kDa band observed in nuclear lysates of cell line HeLa (calculated MW of 46.9kDa according to NP_005216.1). Recommended concentration: 1-3µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES



