

GOAT ANTI-TRANSMEMBRANE PROTEIN 37 ANTIBODY

SKU: EB08033



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	voltage-dependent calcium channel gamma subunit-like protein PR1 PR CACNG5 transmembrane protein 37 TMEM37
Accession ID	NP_899063.2
Blocking Peptide	EBP08033
Immunogen	Peptide with sequence TNQTICFRDLGQ, from the internal region of the protein sequence according to NP_899063.2.
Peptide Sequence	TNQTICFRDLGQ
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	140738
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:32000.
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)