



Telephone: (650) 697-3600

## **GOAT ANTI-TRANSMEMBRANE PROTEIN 37 ANTIBODY**

**SKU:** EB08033

## **SPECIFICATIONS**

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

**Unit Size** 

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

voltage-dependent calcium channel gamma subunit-like protein|PR1|PR|CACNG5|transmembrane protein **Alias** 

37|TMEM37 **Names** 

Accession NP 899063.2 ID

**Blocking** 

EBP08033 **Peptide** 

Peptide with sequence TNQTICFRDLGQ, from the internal region of the protein sequence according to Immunogen

NP 899063.2.

**Peptide** Sequence

**TNQTICFRDLGQ** 

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring medium.png

Method the immunizing peptide.

**Shipping** 

Refrigerated

Instructions

**Predicted Species** 

Human

Human

140738 Gene ID

**Product** 

Grade **ELISA** 

**Detection** Antibody detection limit dilution 1:32000.

Limit

**Application** Pep-ELISA **Type**