


GOAT ANTI-TRC8 ANTIBODY

SKU: EB05902



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym / Alias Names multiple membrane spanning receptor TRC8|patched related protein translocated in renal cancer|RCA1|MGC31961|HRCA1|TRC8|ring finger protein 139|RNFI39

Accession ID NP_009149.2

Blocking Peptide EBP05902

Immunogen Peptide with sequence C-HTGAAAEFNDT, from the C Terminus of the protein sequence according to NP_009149.2.

Peptide Sequence C-HTGAAAEFNDT

Purification Method Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Shipping Instructions Refrigerated

Predicted Species Human, Dog

Reactive Species Human

Human Gene ID 11236

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

IHC Results In paraffin embedded Human Placenta shows membranous staining in the cytotrophoblasts. Recommended concentration, 3-4µg/ml.

ELISA Detection Limit Antibody detection limit dilution 1:16000.

Western Blot Approx 75kDa band observed in Human Skeletal Muscle lysates (calculated MW of 76.0kDa according to NP_009149.2). Recommended concentration: 0.3-1µg/ml.

Application Type Pep-ELISA, WB, IHC

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

