



Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **GOAT ANTI-TRC8 ANTIBODY**

**SKU:** EB05902



Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600


250kDa 150kDa 100kDa 75kDa

50kDa

37kDa

25kDa

20kDa





Telephone: (650) 697-3600

## **SPECIFICATIONS**

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

**Unit Size** 100 μg

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias

multiple membrane spanning receptor TRC8|patched related protein translocated in renal

cancer|RCA1|MGC31961|HRCA1|TRC8|ring finger protein 139|RNF139

**Accession** 

Names

NP 009149.2

**Blocking** 

**Peptide** 

EBP05902

Peptide with sequence C-HTGAAAEEFNDDT, from the C Terminus of the protein sequence according to **Immunogen** 

NP 009149.2.

**Peptide** 

Sequence

C-HTGAAAEEFNDDT

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

**Shipping** 

Instructions

Refrigerated

**Predicted** 

**Species** 

Human, Dog

Reactive

**Species** 

Human

Human

11236 Gene ID

**Product** 

Grade

 $https://prod-vector-labs-pimcore-assets.s 3. us-east-1. a mazon aws.com/assets/products/image/elite\_medium.png$ 

**IHC Results** 

In paraffin embedded Human Placenta shows membranous staining in the cytotrophoblasts. Recommended

concentration, 3-4µg/ml.

**ELISA** 

**Detection** 

Antibody detection limit dilution 1:16000.

Limit

Western Approx 75kDa band observed in Human Skeletal Muscle lysates (calculated MW of 76.0kDa according to

**Blot** NP\_009149.2). Recommended concentration: 0.3-1µg/ml.

**Application** 

Type

Pep-ELISA, WB, IHC

## **GALLERY IMAGES**



Telephone: (650) 697-3600



250kDa 150kDa 100kDa 75kDa

50kDa

37kDa

25kDa

20kDa

