

GOAT ANTI-TRF4-2 / PAPD5 ANTIBODY

SKU: EB05590



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias	FLJ40270 PAP associated domain containing 5 topoisomerase-related function protein 4-2 TRF4-2 PAPD5
Names	
Accession ID	BAA91518.1
Blocking Peptide	EBP05590
Immunogen	Peptide with sequence QCRPLLPHSPALVT, from the C Terminus of the protein sequence according to BAA91518.1.
Product Comments	The NCBI RefSeq NP_071892 has been withdrawn. The immunizing peptide still corresponds to protein BAA91518.1 in the NCBI database.
Peptide Sequence	QCRPLLPHSPALVT
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	64282
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:1000.
Western Blot	Approx 20kDa band seen in Jurkat lysate [Predicted MW of 20kDa according to NP_071892]. Recommended for use at 1-3µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

