


GOAT ANTI-TRIM11 (AA249-263) ANTIBODY

SKU: EB11287



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym / Alias Names TRIM11| tripartite motif-containing protein 11| tripartite motif-containing 11| tripartite motif containing 11| RNF92| RING finger protein 92| E3 ubiquitin-protein ligase TRIM11| BIA1

Accession ID NP_660215.1

Blocking Peptide EBP11287

Immunogen Peptide with sequence C-DALRRVQDVKLPPE, from the internal region of the protein sequence according to NP_660215.1.

Peptide Sequence C-DALRRVQDVKLPPE

Purification Method Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Shipping Instructions Refrigerated

Predicted Species	Human, Pig, Cow
Reactive Species	Human
Human Gene ID	81559
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
IHC Results	Paraffin embedded Human Skeletal Muscle. Recommended concentration: 5µg/ml.
ELISA	
Detection Limit	Antibody detection limit dilution 1:32000.
Western Blot	Approx 52kDa band observed in Human Spleen lysates (calculated MW of 52.8kDa according to NP_660215.1). Recommended concentration: 0.3-1µg/ml. Primary incubation was 1 hour.
Application Type	Pep-ELISA, WB, IHC

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

