


Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-TRIM11 ANTIBODY

SKU: EB05908



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias **Names** tripartite motif-containing 11|RNF92|BIA1|TRIM11

Accession

NP 660215.1

Blocking

EBP05908

Peptide

Immunogen

Peptide with sequence CRPKGGSGDTLAPQ, from the C Terminus of the protein sequence according to NP 660215.1.

Peptide

Sequence

CRPKGGSGDTLAPQ

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.





Telephone: (650) 697-3600



Shipping Instructions

Refrigerated

Predicted

Human

Species Reactive

Human

Species Human

Gene ID

81559

Product

 $https://prod-vector-labs-pimcore-assets.s 3. us-east-1. amazonaws.com/assets/products/image/elite_medium.png$ Grade

In paraffin embedded Human Kidney shows cytoplasm staining preferentially in PCT over DCT. Recommended **IHC Results**

concentration, 4-6µg/ml.

ELISA

Detection

Antibody detection limit dilution 1:32000.

Limit

Application

Type

Pep-ELISA, IHC

GALLERY IMAGES

