

GOAT ANTI-TRIM2 (AA233-247) ANTIBODY

SKU: EB10871

250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym / Alias Names tripartite motif-containing protein 2| tripartite motif protein TRIM2| tripartite motif protein 2| RING finger protein 86| RNF86| KIAA0517| tripartite motif-containing 2|TRIM2

Accession ID NP_056086.2; NP_001123539.1

Blocking Peptide EBP10871

Immunogen Peptide with sequence C-STFDELQKTLNVRKS, from the internal region of the protein sequence according to NP_056086.2; NP_001123539.1.

Product Comments This antibody is expected to recognize both reported isoforms (NP_056086.2; NP_001123539.1).

Peptide Sequence C-STFDELQKTLNVRKS

Purification Method Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Shipping Instructions Refrigerated

Predicted Species Human, Mouse, Rat, Dog, Cow

Reactive Species Mouse

Human Gene ID 23321

Mouse Gene ID 80890

Rat Gene ID 361970

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

ELISA Detection Limit Antibody detection limit dilution 1:16000.

Western Blot Approx 80kDa band observed in Mouse Brain lysates (calculated MW of 81.4kDa according to Mouse NP_109631.1). Recommended concentration: 0.1-0.3µg/ml.

Application Type Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

