

GOAT ANTI-TRIM36 / RNF98 ANTIBODY

SKU: EB05930



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	TRIM36 RNF98 RBCC728 tripartite motif-containing 36 tripartite motif protein 36 zinc-binding protein Rbcc728 HAPRIN RING finger protein 98
Accession ID	NP_061170.2
Blocking Peptide	EBP05930
Immunogen	Peptide with sequence C-PITAKYLEYQEDM, from the C Terminus of the protein sequence according to NP_061170.2.
Peptide Sequence	C-PITAKYLEYQEDM
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat
Reactive Species	Human
Human Gene ID	55521
Mouse Gene ID	28105
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:16000.
Western Blot	Approx 85kDa band observed in Human Brain (Amygdala) lysates (calculated MW of 83.7kDa according to NP_061170.2). Recommended concentration: 1-3µg/ml. Primary incubation was 1 hour.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

