

GOAT ANTI-TRIM6 ANTIBODY

SKU: EB05936



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias Names	RNF89 tripartite motif-containing 6 TRIM6
Accession ID	NP_001003818.1; NP_477514.1
Blocking Peptide	EBP05936
Immunogen	Peptide with sequence CVIPMTLRRPSS, from the C Terminus of the protein sequence according to NP_001003818.1; NP_477514.1.
Product Comments	This antibody is expected to recognise both reported isoforms.
Peptide Sequence	CVIPMTLRRPSS
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	117854
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:2000.
Western Blot	Two bands observed between 25-30kDa in Human Liver and Kidney lysates at 1ug of primary. Please note that the band we observe is different from the predicted size of approx. 56kDa according to NP_477514. Currently we cannot find an explanation for this in the literature. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates? Have any splice variants/modified forms been reported?
Application Type	Pep-ELISA