

GOAT ANTI-TRIP4 / ASC1 (AA94-106) ANTIBODY

SKU: EB12060

SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	thyroid receptor-interacting protein 4 thyroid receptor interacting protein 4 activating signal cointegrator 1 TRIP-4 TR-interacting protein 4 HsT17391 ASC1 ASC-1 thyroid hormone receptor interactor 4 TRIP4
Accession ID	NP_057297.2
Blocking Peptide	EBP12060
Immunogen	Peptide with sequence C-DGQKSGDHLKRGR, from the internal region of the protein sequence according to NP_057297.2.
Peptide Sequence	C-DGQKSGDHLKRGR
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	9325
Mouse Gene ID	56404
Rat Gene ID	315769
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:128000.
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)