Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

## **GOAT ANTI-TRIP4 / ASC1 (AA94-106) ANTIBODY**

**SKU:** EB12060

## **SPECIFICATIONS**

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

**Unit Size** 

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

thyroid receptor-interacting protein 4|thyroid receptor interacting protein 4|activating signal cointegrator **Alias** 1|TRIP-4|TR-interacting protein 4|HsT17391|ASC1|ASC-1|thyroid hormone receptor interactor 4|TRIP4 **Names** 

Accession

NP 057297.2 ID

**Blocking** 

EBP12060 **Peptide** 

Peptide with sequence C-DGQKSGDHLKRGR, from the internal region of the protein sequence according to Immunogen

NP 057297.2.

**Peptide** 

C-DGQKSGDHLKRGR Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

**Shipping** 

Refrigerated Instructions

**Predicted Species** 

Human

Human

9325

Gene ID Mouse

**Gene ID** 

56404

**Rat Gene ID** 315769

**Product** Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring\_medium.png

**ELISA** 

**Detection** Antibody detection limit dilution 1:128000.

Limit

**Application** Pep-ELISA

**Type**