

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-TRPC6 (C-TERM) ANTIBODY

SKU: EB08831





Telephone: (650) 697-3600



250kDa 150kDa 100kDa

75kDa

50kDa

37kDa

25kDa 20kDa

15kDa





Telephone: (650) 697-3600

SPECIFICATIONS

Unit Size

100 µg

Synonym / Alias

transient receptor potential channel 6|TRP6|FSGS2|FLJ14863|FLJ11098|transient receptor potential cation

channel, subfamily C, member 6|TRPC6

Names Accession

NP_004612.2

Blocking Peptide

EBP08831

Immunogen

Peptide with sequence C-EKLSMEPNQEETNR, from the C Terminus of the protein sequence according to

NP 004612.2.

Peptide

C-EKLSMEPNQEETNR Sequence

Shipping

Refrigerated Instructions

Predicted

Species

Human, Cow

Reactive

Species

Human

Human

7225 Gene ID

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite medium.png

Grade

ELISA Detection

Limit

Antibody detection limit dilution 1:128000.

Western

Approx. 110kDa band observed in Human Placenta lysates (calculated MW of 106Da according to

Blot

NP_004612.2). Recommended concentration: 0.03-0.1 μ g/ml.

Application

Pep-ELISA, WB **Type**

GALLERY IMAGES



Telephone: (650) 697-3600



250kDa 150kDa 100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa