

GOAT ANTI-TRPC6 (C-TERM) ANTIBODY

SKU: EB08831



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias Names	transient receptor potential channel 6 TRP6 FSGS2 FLJ14863 FLJ11098 transient receptor potential cation channel, subfamily C, member 6 TRPC6
Accession ID	NP_004612.2
Blocking Peptide	EBP08831
Immunogen	Peptide with sequence C-EKLSMEPNQEETNR, from the C Terminus of the protein sequence according to NP_004612.2.
Peptide Sequence	C-EKLSMEPNQEETNR
Shipping Instructions	Refrigerated
Predicted Species	Human, Cow
Reactive Species	Human
Human Gene ID	7225
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:128000.
Western Blot	Approx. 110kDa band observed in Human Placenta lysates (calculated MW of 106Da according to NP_004612.2). Recommended concentration: 0.03-0.1µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

