

GOAT ANTI-TRPV3 (AA610-622) ANTIBODY

SKU: EB08468



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias Names	vanilloid receptor-related osmotically activated channel protein vanilloid receptor 3 VRL3 transient receptor potential cation channel, subfamily V, member 3 TRPV3
Accession ID	NP_659505.1
Blocking Peptide	EBP08468
Immunogen	Peptide with sequence EKCPKDNKDCSSY, from the internal region of the protein sequence according to NP_659505.1.
Peptide Sequence	EKCPKDNKDCSSY
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	162514
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:128000.
Western Blot	Preliminary experiments in Human Brain (Cerebellum), Human Skin and Rat Spinal Cord lysates gave no specific signal but low background (at antibody concentration up to 1µg/ml). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)

