

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-TRPV3 (AA610-622) ANTIBODY

SKU: EB08468

SPECIFICATIONS

Unit Size

100 μg

Synonym / Alias

vanilloid receptor-related osmotically activated channel protein|vanilloid receptor 3|VRL3|transient receptor

potential cation channel, subfamily V, member 3|TRPV3

Names Accession

ID

NP_659505.1

Blocking Peptide

EBP08468

Peptide with sequence EKCPKDNKDCSSY, from the internal region of the protein sequence according to **Immunogen**

NP 659505.1.

Peptide

Sequence

EKCPKDNKDCSSY

Shipping

Refrigerated Instructions

Predicted

Species

Human

Human

Gene ID

162514

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

Grade **ELISA**

Detection

Antibody detection limit dilution 1:128000.

Limit

Western Blot

Preliminary experiments in Human Brain (Cerebellum), Human Skin and Rat Spinal Cord lysates gave no specific signal but low background (at antibody concentration up to 1µg/ml). We would appreciate any feedback from

people in the field - have any results been reported with other antibodies/lysates?

Application

Type

Pep-ELISA