

GOAT ANTI-TSNAX ANTIBODY

SKU: EB07118

SPECIFICATIONS

Unit Size	100 µg
Synonym /	
Alias	translin-like protein TRAX translin-associated factor X TSNAX
Names	
Accession ID	NP_005990.1
Blocking Peptide	EBP07118
Immunogen	Peptide with sequence C-KENKTPSSDAQDKQ, from the internal region of the protein sequence according to NP_005990.1.
Peptide Sequence	C-KENKTPSSDAQDKQ
Shipping Instructions	Refrigerated
Predicted Species	Human, Dog
Human Gene ID	7257
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA	
Detection Limit	Antibody detection limit dilution 1:1000.
Western Blot	Preliminary experiments in Human Brain gave no specific signal but low background (at antibody concentration up to 1µg/ml). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?
Application Type	Pep-ELISA