

GOAT ANTI-TSPYL2 (AA2-13) ANTIBODY

SKU: EB11280



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	TSPYL2 TSPY-like protein 2 TSPY-like 2 TSPX testis-specific Y-encoded-like protein 2 testis-specific protein Y encoded-like 2 SE204 nuclear protein of 79 kDa NP79 HRIHFB2216 differentially-expressed nucleolar TGF-beta1 target protein differentially expressed nucleolar TGF-beta1 target DENTT cutaneous T-cell lymphoma-associated tumor antigen se20-4 cutaneous T-cell lymphoma-associated antigen se20-4 CTCL-associated antigen se20-4 CTCL tumor antigen se20-4 CTCL CINAP cell division autoantigen 1 CDA1 CASK-interacting nucleosome assembly protein
Accession ID	NP_071400.1
Blocking Peptide	EBP11280
Immunogen	Peptide with sequence DRPDEGPPAKTR-C, from the N Terminus of the protein sequence according to NP_071400.1.
Peptide Sequence	DRPDEGPPAKTR-C
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Pig, Cow
Human Gene ID	64061
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:4000.
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)