

GOAT ANTI-TSPYL2 (AA441-454) ANTIBODY

SKU: EB11281



250kDa
150kDa
100kDa
75kDa
50kDa
37kDa
25kDa
20kDa

SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	TSPYL2 TSPY-like protein 2 TSPY-like 2 TSPX testis-specific Y-encoded-like protein 2 testis-specific protein Y encoded-like 2 SE204 nuclear protein of 79 kDa NP79 HRIHFB2216 differentially-expressed nucleolar TGF-beta1 target protein differentially expressed nucleolar TGF-beta1 target DENTT cutaneous T-cell lymphoma-associated tumor antigen se20-4 cutaneous T-cell lymphoma-associated antigen se20-4 CTCL-associated antigen se20-4 CTCL tumor antigen se20-4 CTCL CINAP cell division autoantigen 1 CDA1 CASK-interacting nucleosome assembly protein
Accession ID	NP_071400.1
Blocking Peptide	EBP11281
Immunogen	Peptide with sequence C-SDIDETIHDIKISD, from the internal region of the protein sequence according to NP_071400.1.
Peptide Sequence	C-SDIDETIHDIKISD
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	64061
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
IHC Results	Paraffin embedded Human Pancreas. Recommended concentration: 5µg/ml.
ELISA	
Detection Limit	Antibody detection limit dilution 1:64000.
Western Blot	Approx 75kDa band observed in lysates of cell lines Jurkat and MOLT4 (calculated MW of 79.4kDa according to NP_071400.1). Recommended concentration: 0.3-1µg/ml. Primary incubation was 1 hour.
Application Type	Pep-ELISA, WB, IHC

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

