


Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **GOAT ANTI-TTLL4 ANTIBODY**

**SKU:** EB09832



## **SPECIFICATIONS**

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

**Unit Size** 

Storage

Aliquot and store at -20°C. Minimize freezing and thawing. Instructions

Synonym / **Alias** 

KIAA0173|OTTHUMP00000164099|tubulin tyrosine ligase-like family, member 4|TTLL4

**Names** 

Accession

NP\_055455.3

ID **Blocking** 

**Peptide** 

EBP09832

Peptide with sequence C-DSQDPTKEIRFTE, from the internal region of the protein sequence according to **Immunogen** 

NP\_055455.3.

**Peptide** Sequence

C-DSQDPTKEIRFTE

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping Instructions

Refrigerated





Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

**Predicted** 

Human, Mouse, Rat, Cow

**Species** Reactive

Human

**Species** Human

9654

**Gene ID** 

Mouse

67534

**Gene ID** 

**Rat Gene ID** 690512

**Product** 

Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_medium.png

**IHC Results** 

In paraffin embedded Human Lung shows cytoplasm staining in the resporatory epithelium of the bronchus.

Recommended concentration, 5-10µg/ml.

**ELISA** 

Detection

Antibody detection limit dilution 1:32000.

Limit

**Type** 

**Application** Pep-ELISA, IHC

## **GALLERY IMAGES**

