


GOAT ANTI-TTLL4 ANTIBODY

SKU: EB09832



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias KIAA0173|OTTHUMP00000164099|tubulin tyrosine ligase-like family, member 4|TTLL4

Names

Accession ID NP_055455.3

Blocking Peptide EBP09832

Immunogen Peptide with sequence C-DSQDPTKEIRFTE, from the internal region of the protein sequence according to NP_055455.3.

Peptide Sequence C-DSQDPTKEIRFTE

Purification Method Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Shipping Instructions Refrigerated

Predicted Species	Human, Mouse, Rat, Cow
Reactive Species	Human
Human Gene ID	9654
Mouse Gene ID	67534
Rat Gene ID	690512
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
IHC Results	In paraffin embedded Human Lung shows cytoplasm staining in the respiratory epithelium of the bronchus. Recommended concentration, 5-10µg/ml.
ELISA	
Detection Limit	Antibody detection limit dilution 1:32000.
Application Type	Pep-ELISA, IHC

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

