


# GOAT ANTI-TUSC4 / NPRL2 (AA 269 TO 280) ANTIBODY

**SKU:** EB08276



250kDa  
150kDa  
100kDa  
75kDa  
50kDa  
37kDa  
25kDa  
20kDa  
15kDa

---

## SPECIFICATIONS

|                              |                                                                                                                                                                                                                                   |
|------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Unit Size</b>             | 100 µg                                                                                                                                                                                                                            |
| <b>Synonym / Alias Names</b> | homologous to yeast nitrogen permease (candidate tumor suppressor) 2810446G01Rik NPR2L tumor suppressor candidate 4 NPRL2 TUSC4                                                                                                   |
| <b>Accession ID</b>          | NP_006536.3                                                                                                                                                                                                                       |
| <b>Blocking Peptide</b>      | EBP08276                                                                                                                                                                                                                          |
| <b>Immunogen</b>             | Peptide with sequence C-TKQGHKRASLRD, from the internal region of the protein sequence according to NP_006536.3.                                                                                                                  |
| <b>Peptide Sequence</b>      | C-TKQGHKRASLRD                                                                                                                                                                                                                    |
| <b>Shipping Instructions</b> | Refrigerated                                                                                                                                                                                                                      |
| <b>Predicted Species</b>     | Human, Mouse                                                                                                                                                                                                                      |
| <b>Reactive Species</b>      | Human                                                                                                                                                                                                                             |
| <b>Human Gene ID</b>         | 10641                                                                                                                                                                                                                             |
| <b>Mouse Gene ID</b>         | 56032                                                                                                                                                                                                                             |
| <b>Rat Gene ID</b>           | 363138                                                                                                                                                                                                                            |
| <b>Product Grade</b>         | <a href="https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png">https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png</a> |
| <b>ELISA Detection Limit</b> | Antibody detection limit dilution 1:16000.                                                                                                                                                                                        |
| <b>Western Blot</b>          | Approx 48kDa band observed in Human Skeletal Muscle lysates (calculated MW of 43.7kDa according to NP_006536.3). Recommended concentration: 1-3µg/ml.                                                                             |
| <b>Application Type</b>      | Pep-ELISA, WB                                                                                                                                                                                                                     |

## DOCUMENTS

- [Data Sheet](#)

## GALLERY IMAGES

