



www.everestbiotech.com

Email: customerservice@vectorlabs.com

Telephone: [\(650\) 697-3600](tel:(650)697-3600)

GOAT ANTI-TXL2 / TXNDC6 ANTIBODY

SKU: EB08520



250kDa
150kDa
100kDa
75kDa

50kDa
37kDa

25kDa
20kDa

15kDa



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	thioredoxin-like 2 TXL-2 MGC129586 thioredoxin domain containing 6 TXL2 TXNDC6
Accession ID	NP_835231.1
Blocking Peptide	EBP08520
Immunogen	Peptide with sequence C-ESSTQPRLKITDLD, from the C Terminus of the protein sequence according to NP_835231.1.
Peptide Sequence	C-ESSTQPRLKITDLD
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	347736
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:8000.
Western Blot	Approx 38kDa band observed in Human Uterus lysates (calculated MW of 29.7kDa according to NP_835231.1). The observed molecular weight corresponds to earlier findings in literature with different antibodies (Sadek et al, J Biol Chem. 2003 Apr 11;278(15):13133-42; 12569107). Recommended concentration: 1-3µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)



GALLERY IMAGES

