

GOAT ANTI-TXNDC1 / TMX ANTIBODY

SKU: EB05086



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	transmembrane Trx-related protein protein disulfide isomerase family A, member 11 PDIA11 thioredoxin-related transmembrane protein 1 DKFZP564E1962 TXNDC1 thioredoxin domain containing 1 TMX TXNDC TMX1
Accession ID	NP_110382.3
Blocking Peptide	EBP05086
Immunogen	Peptide with sequence C-RSLGPSLATDKS, from the C Terminus of the protein sequence according to NP_110382.3.
Peptide Sequence	C-RSLGPSLATDKS
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	81542
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA	
Detection Limit	Antibody detection limit dilution 1:64000.
Western Blot	Approx 35kDa band observed in Human Liver lysates (calculated MW of 31.8kDa according to NP_110382.3). Recommended concentration: 0.01-0.03µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

