

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-TXNDC3 ANTIBODY

SKU: EB08537

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Synonym /

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

thioredoxin domain-containing 3 (spermatozoa)|thioredoxin domain containing 3|sperm-specific thioredoxin **Alias** 2|NM23-H8|SPTRX2|NME8|CILD6|thioredoxin domain containing 3 (spermatozoa)|TXNDC3

Names

Accession

NP 057700.3 ID

Blocking

EBP08537 **Peptide**

Peptide with sequence C-TEPQTDTEPNERSE, from the internal region of the protein sequence according to Immunogen

NP 057700.3.

Peptide

C-TEPQTDTEPNERSE Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Species

Human

Human

51314 Gene ID

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring medium.png

Grade **ELISA**

Detection

Antibody detection limit dilution 1:16000.

Limit

Application Pep-ELISA **Type**