

GOAT ANTI-TXNDC3 ANTIBODY

SKU: EB08537



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	thioredoxin domain-containing 3 (spermatozoa) thioredoxin domain containing 3 sperm-specific thioredoxin 2 NM23-H8 SPTRX2 NME8 CILD6 thioredoxin domain containing 3 (spermatozoa) TXNDC3
Accession ID	NP_057700.3
Blocking Peptide	EBP08537
Immunogen	Peptide with sequence C-TEPQTDTEPNERSE, from the internal region of the protein sequence according to NP_057700.3.
Peptide Sequence	C-TEPQTDTEPNERSE
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	51314
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:16000.
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)