



www.everestbiotech.com

Email: customerservice@vectorlabs.com

Telephone: [\(650\) 697-3600](tel:(650)697-3600)

GOAT ANTI-U-PAR / CD87 ANTIBODY

SKU: EB12639





SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	PLAUR plasminogen activator, urokinase receptor CD87 U-PAR UPAR URKR monocyte activation antigen Mo3 u-plasminogen activator receptor form 2 urokinase plasminogen activator surface receptor urokinase-type plasminogen activator (uPA) receptor
Accession ID	NP_002650.1; NP_001005376.1; NP_001005377.1
Blocking Peptide	EBP12639
Immunogen	Peptide with sequence C-RTGLKITSLTE , from the internal region of the protein sequence according to NP_002650.1; NP_001005376.1; NP_001005377.1.
Product Comments	This antibody is expected to recognize all reported isoforms (NP_002650.1; NP_001005376.1; NP_001005377.1).
Peptide Sequence	C-RTGLKITSLTE
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	5329
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:128000.
Western Blot	Approx 60kDa band observed in lysates of cell line U251 (calculated MW of 37.0kDa according to NP_002650.1). The observed molecular weight corresponds to the glycosylated form. Recommended concentration: 1-3µg/ml.
Application Type	Pep-ELISA, WB



www.everestbiotech.com

Email: customerservice@vectorlabs.com

Telephone: [\(650\) 697-3600](tel:(650)697-3600)

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

