

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **GOAT ANTI-U-PAR / CD87 ANTIBODY**

**SKU:** EB12639



Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

250kDa 150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa





Telephone: (650) 697-3600

## **SPECIFICATIONS**

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

**Unit Size** 100 μg

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

PLAUR|plasminogen activator, urokinase receptor|CD87|U-PAR|UPAR|URKR|monocyte activation antigen Mo3|uplasminogen activator receptor form 2|urokinase plasminogen activator surface receptor|urokinase-type

**Alias** Names

plasminogen activator (uPA) receptor

**Accession** 

NP\_002650.1; NP\_001005376.1; NP\_001005377.1

**Blocking** 

**Peptide** 

EBP12639

**Immunogen** 

Peptide with sequence C-RTGLKITSLTE, from the internal region of the protein sequence according to

NP\_002650.1; NP\_001005376.1; NP\_001005377.1.

**Product** This antibody is expected to recognize all reported isoforms (NP\_002650.1; NP\_001005376.1;

NP\_001005377.1). Comments

**Peptide** 

Sequence

C-RTGLKITSLTE

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Human

Human

5329

Shipping

Refrigerated Instructions

**Predicted** 

**Species** 

Reactive **Species** 

Human

Gene ID

**Product** 

Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite medium.png

**ELISA** 

**Detection** 

Antibody detection limit dilution 1:128000.

Limit

Blot

Western

Approx 60kDa band observed in lysates of cell line U251 (calculated MW of 37.0kDa according to NP\_002650.1). The observed molecular weight corresponds to the glycosylated form. Recommended

concentration: 1-3µg/ml.

**Application** 

Pep-ELISA, WB **Type** 

## **GALLERY IMAGES**



Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

250kDa 150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa