

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

## **GOAT ANTI-UBE1L2 ANTIBODY**

**SKU:** EB08151

## **SPECIFICATIONS**

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

**Unit Size** 

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

UBA6, ubiquitin-activating enzyme E1|UBA6|MOP-4|FLJ23367|FLJ10808|ubiquitin-activating enzyme E1-like **Alias** 

2|UBE1L2 **Names** 

Accession

NP 060697.4

Peptide with sequence C-ETTEVRKTKIRN, from the internal region of the protein sequence according to Immunogen NP\_060697.4.

**Peptide** 

C-ETTEVRKTKIRN Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping

Refrigerated Instructions

**Predicted** 

Human **Species** 

Human

55236 Gene ID

**Product** 

Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring\_medium.png

**ELISA** 

**Detection** Antibody detection limit dilution 1:4000.

Limit

Western Blot: Preliminary experiments in Human Testis and Liver lysates gave no specific signal but low

Western background (at antibody concentration up to 1µg/ml). We would appreciate any feedback from people in the field -**Blot** 

have any results been reported with other antibodies/lysates?

Application

Pep-ELISA **Type**