

GOAT ANTI-UBE1L2 ANTIBODY

SKU: EB08151



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	UBA6, ubiquitin-activating enzyme E1 UBA6 MOP-4 FLJ23367 FLJ10808 ubiquitin-activating enzyme E1-like 2 UBE1L2
Accession ID	NP_060697.4
Immunogen	Peptide with sequence C-ETTEVRKTKIRN, from the internal region of the protein sequence according to NP_060697.4.
Peptide Sequence	C-ETTEVRKTKIRN
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	55236
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:4000.
Western Blot	Western Blot: Preliminary experiments in Human Testis and Liver lysates gave no specific signal but low background (at antibody concentration up to 1µg/ml). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?
Application Type	Pep-ELISA