


Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-UBR2 ANTIBODY

SKU: EB09658



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin

Unit Size 100 μg Storage

Storage Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym / ubiquitin ligase E3 alpha-

II|OTTHUMP00000039768|OTTHUMP00000039767|OTTHUMP00000016405|OTTHUMP00000016403|d)392M17.3|d)242G1.1|bA49A4.1|RP3-392M17.3|MGC71112|KIAA0349|DKFZp686C08114|C6orf133|ubiquitin

Names protein ligase E3 component n-recognin 2|UBR2

Accession NP_056070.1 Blocking Peptide

Immunogen Peptide with sequence SCPKSGGDKSRAPT, from the internal region of the protein sequence according to NP_056070.1.

Purification
Mathed

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Shipping Instructions Refrigerated

Predicted Human, Mouse, Dog, Cow

Reactive Human 23304

Product

IHC Results Paraffin embedded Brain (Cortex), Recommended concentration; 5µg/ml. ELISA

Detection Antibody detection limit dilution 1:32000.

Preliminary experiments gave bands at approx 150kDa in Rat Testis lysates after 0.01µg/ml antibody staining. Please note that currently we cannot find an explanation in the literature for the bands we observe given the calculated size of 199kDa according to Rat NP_001171542.1. The detected band was successfully blocked by incubation with the immunizing peptide (and BLAST results with the immunizing peptide Western Blot

sequence did not identify any other proteins to explain the additional bands).





Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

Application Pep-ELISA, IHC

GALLERY IMAGES

