




Telephone: (650) 697-3600



## **GOAT ANTI-UBR4 / ZUBR1 ANTIBODY**

**SKU:** EB09611



## **SPECIFICATIONS**

**Formulation** Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 μc

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

-

**Synonym** / UBR4|ubiquitin protein ligase E3 component n-recognin

Alias

 $4|FLJ41863|KIAA0462|KIAA1307|RBAF600|ZUBR1|p600|retinoblastoma-associated\ factor\ 600|retinoblastoma-associated\ factor\ 600|retinob$ 

Names associated factor 600-like protein|zinc finger, UBR1 type 1

Accession

NP\_065816.2

Blocking

ID

EBP09611

Peptide

Peptide with sequence C-NTTRATSREEKNLQ, from the C Terminus of the protein sequence according to

NP\_065816.2.

Peptide

C-NTTRATSREEKNLQ

Sequence Purification

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

**Method** using the immunizing peptide.

Shipping

Instructions Refrigerated





Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

**Predicted** 

Human, Mouse, Rat, Dog, Cow

Species Reactive

Human

Species Human

Gene ID

23352

Mouse

Gene ID

69116

Rat Gene ID 313658

Product Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_medium.png

**IHC Results** 

In paraffin embedded Spleen shows various types of subnuclear staining in select splenocytes and additional

membrane plus cytoplasm staining elsewhere. Recommended concentration, 5-10µg/ml.

**ELISA** 

Detection Limit Antibody detection limit dilution 1:128000.

Western

Not yet tested - our routinely used western blotting protocol does not allow detection of proteins as large as the calculated size of 574kDa according to NP\_065816.2. Therefore we cannot recommend an optimal concentration and the antibody is an aspiring product. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates? Have any further splice variants/modified

forms been reported?

Application

Type

**Blot** 

Pep-ELISA, IHC

## **GALLERY IMAGES**

