

EBP05415 - IMMUNIZING PEPTIDE

SKU: EBP05415



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	Shipped at ambient temperature, store at -20°C.
Instructions	
Accession ID	NP_068823.2; NP_954595.1
Immunogen	Peptide with sequence PGEVQASYLKSQS-C, from the N Terminus of the protein sequence according to NP_068823.2; NP_954595.1.
Product Comments	The peptide was used in the production of EB05415 - Goat Anti-UEV1 / CROC1A Antibody. It is currently untested in blocking applications.
Peptide Sequence	PGEVQASYLKSQS-C
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
Western Blot	Approx 18-23kDa band seen in Human Brain lysate. Recommended for use at 0.2-1µg/ml.

DOCUMENTS

- [Data Sheet](#)