

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP05415 - IMMUNIZING PEPTIDE

SKU: EBP05415

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

Unit Size 100 µg

Storage

Shipped at ambient temperature, store at -20°C. Instructions

Accession

NP 068823.2; NP 954595.1 ID

Peptide with sequence PGEVQASYLKSQS-C, from the N Terminus of the protein sequence according to Immunogen

NP 068823.2; NP 954595.1.

The peptide was used in the production of EB05415 - Goat Anti-UEV1 / CROC1A Antibody. It is currently **Product**

Comments untested in blocking applications.

Peptide

PGEVQASYLKSQS-C Sequence

Shipping

Instructions Ambient

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

Grade Western

Approx 18-23kDa band seen in Human Brain lysate. Recommended for use at 0.2-1µg/ml. **Blot**