

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-UEV1 / CROC1A ANTIBODY

SKU: EB05415



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa



Telephone: (650) 697-3600

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 μg

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

DNA-binding protein|CROC-1A|UEV1A|UEV-1|UBE2V|CROC1A|CIR1|ubiquitin-conjugating enzyme E2 variant

Alias 1|UEV1|UBE2V1 Names

Accession

NP_068823.2; NP_954595.1

Blocking

EBP05415 **Peptide**

Peptide with sequence PGEVQASYLKSQS-C, from the N Terminus of the protein sequence according to **Immunogen**

NP_068823.2; NP_954595.1.

Product Comments

This antibody is expected to recognize isoform a (NP_068823.2; NP_954595.1) only.

Peptide

Sequence

PGEVQASYLKSQS-C

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Species

Human

Reactive **Species**

Human

Human

7335 Gene ID

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png Grade

ELISA

Detection Antibody detection limit dilution 1:32000.

Limit

Western

Approx 18-23kDa band seen in Human Brain lysate. Recommended for use at 0.2-1µg/ml. **Blot**

Application

Pep-ELISA, WB **Type**

GALLERY IMAGES



Telephone: (650) 697-3600



250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa