

GOAT ANTI-UEV1 / CROC1A ANTIBODY

SKU: EB05415



250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	DNA-binding protein CROC-1A UEV1A UEV-1 UBE2V CROC1A CIR1 ubiquitin-conjugating enzyme E2 variant 1 UEV1 UBE2V1
Accession ID	NP_068823.2; NP_954595.1
Blocking Peptide	EBP05415
Immunogen	Peptide with sequence PGEVQASYLKSQS-C, from the N Terminus of the protein sequence according to NP_068823.2; NP_954595.1.
Product Comments	This antibody is expected to recognize isoform a (NP_068823.2; NP_954595.1) only.
Peptide Sequence	PGEVQASYLKSQS-C
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	7335
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:32000.
Western Blot	Approx 18-23kDa band seen in Human Brain lysate. Recommended for use at 0.2-1µg/ml.
Application Type	Pep-ELISA, WB

GALLERY IMAGES

