

GOAT ANTI-UEV2V2 / MMS2 ANTIBODY

SKU: EB05416



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	UBE2V2 MMS2 ubiquitin-conjugating enzyme E2 variant 2 UEV2 EDPF1 UEV-2 DDVIT1 EDAF-1 EDPF-1 DDVit-1 1 alpha, 25-dihydroxyvitamin D3-inducible enterocyte differentiation promoting factor methyl methanesulfonate sensitive 2, S. cerevisiae, homolog of ubiquitin-conjugating enzyme E2v2
Accession ID	NP_003341.1; NP_954673.1
Blocking Peptide	EBP05416
Immunogen	Peptide with sequence AVSTGVKVP RNF-C, from the N Terminus of the protein sequence according to NP_003341.1; NP_954673.1.
Product Comments	This antibody is also expected to recognize the readthrough transcript product KUA-UEV (NP_954673.1, GeneID 387522). Also cross-reactivity is expected with UEV2V1 (GeneID 7335).
Peptide Sequence	AVSTGVKVP RNF-C
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Dog, Cow
Reactive Species	Human
Human Gene ID	7336
Mouse Gene ID	70620
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:32000.
Western Blot	Approx 20kDa band seen in human brain lysate. Recommended for use at 0.2-1µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

