

GOAT ANTI-UHMK1 ANTIBODY

SKU: EB09397



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	kinase interacting with leukemia-associated gene (stathmin) kinase interacting stathmin OTTHUMP00000029529 KIS protein kinase Kist KIS U2AF homology motif (UHM) kinase 1 UHMK1
Accession ID	NP_787062.1
Blocking Peptide	EBP09397
Immunogen	Peptide with sequence C-AGDSKAAQKLLTGR, from the C-Terminus (near) of the protein sequence according to NP_787062.1.
Peptide Sequence	C-AGDSKAAQKLLTGR
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Dog, Cow
Human Gene ID	127933
Mouse Gene ID	16589
Rat Gene ID	246332
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:2000.
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)