Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-UHMK1 ANTIBODY

SKU: EB09397

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

kinase interacting with leukemia-associated gene (stathmin)|kinase interacting **Alias**

stathmin|OTTHUMP00000029529|KIS protein kinase|Kist|KIS|U2AF homology motif (UHM) kinase 1|UHMK1

Names

Accession NP 787062.1

Blocking

EBP09397 **Peptide**

Peptide with sequence C-AGDSKAAQKLLTGR, from the C-Terminus (near) of the protein sequence according to Immunogen

NP 787062.1.

Peptide

C-AGDSKAAQKLLTGR Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Species

Human, Mouse, Rat, Dog, Cow

Human Gene ID

127933

Mouse **Gene ID**

16589

Rat Gene ID 246332

Product Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

ELISA

Detection Antibody detection limit dilution 1:2000.

Limit

Application Pep-ELISA

Type