

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-UPF3A / RENT3A ANTIBODY

SKU: EB09661

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

UPF3 regulator of nonsense transcripts homolog A|OTTHUMP00000018789|UPF3|RENT3A|HUPF3A|UPF3 regulator **Alias**

of nonsense transcripts homolog A (yeast)|UPF3A **Names**

Accession

NP 075387.1; NP 542418.1 ID

Blocking

EBP09661 **Peptide**

Peptide with sequence C-KETDKQKKIAEKE, from the internal region of the protein sequence according to Immunogen

NP 075387.1; NP 542418.1.

Product

This antibody is expected to recognize both isoforms (NP_075387.1; NP_542418.1). **Comments**

Peptide

Sequence

C-KETDKQKKIAEKE

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted Species

Human, Cow

Human

65110

Product

Gene ID

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade

ELISA

Detection Antibody detection limit dilution 1:8000.

Limit

Western

Preliminary experiments gave an approx 48kDa band in lysates of cell lines HeLa and Kelly after 0.3µg/ml antibody staining. Please note that currently we cannot find an explanation in the literature for the band we observe given

the calculated size of 51kDa according to NP 542418.1 and the reported sizes of 65 to 70kDa from different

sources. The 48kDa band was successfully blocked by incubation with the immunizing peptide.

Application

Type

Blot

Pep-ELISA





Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600