

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600


GOAT ANTI-UROCORTIN 3 ANTIBODY

SKU: EB08308



Telephone: (650) 697-3600





Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage

Aliquot and store at -20°C. Minimize freezing and thawing. Instructions

Synonym /

Alias urocortin III|urocortin 3|stresscopin|UCNIII|SPC|SCP|MGC119002|urocortin 3 (stresscopin)|UCN3

Names

Accession NP 444277.2 ID

Blocking

EBP08308 **Peptide**

Peptide with sequence C-KPRQDTAKSPHRTK, from the internal region of the protein sequence according to **Immunogen**

NP 444277.2.

Peptide

C-KPRQDTAKSPHRTK Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Human **Species**

Reactive

Human

Species Human

114131

Gene ID Product

Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

ELISA

Detection Antibody detection limit dilution 1:16000.

Limit

Western **Blot**

In transfected HEK293 transiently expressing full-length Human UCN3 (myc and DYKDDDDK tagged), a band of approx. 19kDa was observed. No bands were observed in mock-transfected HEK293 and the same band was

observed using anti-DYKDDDDK tag antibody. Recommended concentration, 1- 2μg/ml.

Application

Type

Pep-ELISA, WB

GALLERY IMAGES





Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

