


GOAT ANTI-UROCORTIN 3 ANTIBODY

SKU: EB08308



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	urocortin III urocortin 3 stresscopin UCNIII SPC SCP MGC119002 urocortin 3 (stresscopin) UCN3
Accession ID	NP_444277.2
Blocking Peptide	EBP08308
Immunogen	Peptide with sequence C-KPRQDTAKSPHRTK, from the internal region of the protein sequence according to NP_444277.2.
Peptide Sequence	C-KPRQDTAKSPHRTK
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	114131
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:16000.
Western Blot	In transfected HEK293 transiently expressing full-length Human UCN3 (myc and DYKDDDDK tagged), a band of approx. 19kDa was observed. No bands were observed in mock-transfected HEK293 and the same band was observed using anti-DYKDDDDK tag antibody. Recommended concentration, 1- 2µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

