


GOAT ANTI-UROMODULIN (AA525-536) ANTIBODY

SKU: EB11647



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym / Alias uromucoid|uromodulin (uromucoid, Tamm-Horsfall glycoprotein)|uromodulin|UMOD|THP|THGP|tamm-Horsfall urinary glycoprotein|Tamm-Horsfall

Names glycoprotein|OTTHUMP00000162213|OTTHUMP00000162212|MCKD2|HNFJ1|HNFJ|FJHN|ADMCKD2

Usage Immunofluorescence: Strong expression of the protein seen in the Plasma membrane and

Summary Golgi of HeLa cells. Recommended concentration: 10µg/ml.

Accession ID NP_003352.2

Blocking Peptide EBP11647

Immunogen Peptide with sequence DRCPHTRDSTIQ, from the internal region of the protein sequence according to NP_003352.2.

Product Comments	Reported variants represent identical protein: NP_001008390.1, NP_003352.2
Peptide Sequence	DRCPHTRDSTIQ
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	7369
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:64000.
Application Type	Pep-ELISA, IF

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

